


HSP10 Antibody (R34993)

Catalog No.	Formulation	Size
R34993-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3336
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HSP10 antibody is available for research use only.



Western blot testing of HepG2 lysate with HSP10 antibody at 0.01 ug/ml. Expected molecular weight: ~10kDa.

Description

Additional name(s) for this target protein: Heat shock protein 10; HSPE1

Application Notes

Optimal dilution of the HSP10 antibody should be determined by the researcher.

Immunogen

Amino acids PEKSQGKVLQAT were used as the immunogen for this HSP10 antibody.

Storage

Aliquot and store the HSP10 antibody at -20°C.