

## HSP22 Antibody (R34788)

Catalog No.	Formulation	Size
R34788-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26353
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HSP22 antibody is available for research use only.



Western blot testing of human muscle lysate with HSP22 antibody at 1ug/ml. Observed molecular weight: 22-27 kDa.

## Description

Additional name(s) for this target protein: Heat shock 22kDa protein 8; Heat shock 27kDa protein 8; H11; HSPB8

## Application Notes

Optimal dilution of the HSP22 antibody should be determined by the researcher.

## Immunogen

Amino acids NELPQDSQEVTCT were used as the immunogen for this HSP22 antibody.

## Storage

Aliquot and store the HSP22 antibody at -20oC.