


HSP22 Antibody (R34788)

Catalog No.	Formulation	Size
R34788-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26353
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HSP22 antibody is available for research use only.



Western blot testing of human muscle lysate with HSP22 antibody at 1ug/ml. Observed molecular weight: 22-27 kDa.

Description

Additional name(s) for this target protein: Heat shock 22kDa protein 8; Heat shock 27kDa protein 8; H11; HSPB8

Application Notes

Optimal dilution of the HSP22 antibody should be determined by the researcher.

Immunogen

Amino acids NELPQDSQEVCT were used as the immunogen for this HSP22 antibody.

Storage

Aliquot and store the HSP22 antibody at -20°C.