


HSPA13 Antibody (R34006)

Catalog No.	Formulation	Size
R34006-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6782
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HSPA13 antibody is available for research use only.



Western blot testing of human uterus lysate with HSPA13 antibody at 0.1ug/ml.
Predicted molecular weight: ~52kDa, routinely observed at ~70kDa.

Description

Additional name(s) for this target protein: Heat shock protein 70kDa family, member 13, STCH

Application Notes

Optimal dilution of the HSPA13 antibody should be determined by the researcher.

Immunogen

Amino acids ESVELADSNPQNT were used as the immunogen for this HSPA13 antibody.

Storage

Aliquot and store the HSPA13 antibody at -20oC.