


HSPC117 Antibody (R35184)


Catalog No.	Formulation	Size
R35184-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51493
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This HSPC117 antibody is available for research use only.



Western blot testing of HeLa lysate with HSPC117 antibody at 0.3ug/ml.
Predicted/observed molecular weight: ~55kDa.



IHC staining of FFPE human testis with HSPC117 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: C22orf28

Application Notes

Optimal dilution of the HSPC117 antibody should be determined by the researcher.

Immunogen

Amino acids QADPNKVSARAKKR were used as the immunogen for this HSPC117 antibody.

Storage

Aliquot and store the HSPC117 antibody at -20oC.