


HSPC117 Antibody (R35675)

Catalog No.	Formulation	Size
R35675-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51493
Applications	Western Blot : 0.05-0.2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This HSPC117 antibody is available for research use only.



Western blot testing of Jurkat lysate with HSPC117 antibody at 0.05ug/ml.
Predicted/observed molecular weight: ~55kDa.

Description

Additional name(s) for this target protein: Chromosome 22 open reading frame 28; C22orf28

Application Notes

Optimal dilution of the HSPC117 antibody should be determined by the researcher.

Immunogen

Amino acids RLLRTNLDESVDVQP were used as the immunogen for this HSPC117 antibody.

Storage

Aliquot and store the HSPC117 antibody at -20oC.