


Htr1a Antibody (R33963)

Catalog No.	Formulation	Size
R33963-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3350
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Htr1a antibody is available for research use only.



Western blot testing of mouse brain lysate with Htr1a antibody at 1ug/ml. Predicted molecular weight: ~46 kDa.

Description

Additional name(s) for this target protein: 5-hydroxytryptamine (serotonin) receptor 1A

Application Notes

Optimal dilution of the Htr1a antibody should be determined by the researcher.

Immunogen

Amino acids ERKNERTAEAKRK were used as the immunogen for this Htr1a antibody.

Storage

Aliquot and store the Htr1a antibody at -20oC.