


## HTR1B Antibody / 5HT1B (R35985)

Catalog No.	Formulation	Size
R35985-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cat, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3351
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This HTR1B antibody is available for research use only.



Western blot testing of human brain lysate with HTR1B antibody at 2ug/ml. Predicted molecular weight: ~44 kDa.

## Description

Additional name(s) for this target protein: 5-hydroxytryptamine (serotonin) receptor 1B

## Application Notes

Optimal dilution of the HTR1B antibody should be determined by the researcher.

## Immunogen

Amino acids SNEDFKQAFHK were used as the immunogen for this HTR1B antibody.

## **Storage**

Aliquot and store the HTR1B antibody at -20oC.