


## Htr3b Antibody (R35234)

Catalog No.	Formulation	Size
R35234-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	58963
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Htr3b antibody is available for research use only.



Western blot testing of rat brain lysate with Htr3b antibody at 0.1ug/ml. Predicted molecular weight (rat): ~50kDa.

## Description

Additional name(s) for this target protein: 5-hydroxytryptamine (serotonin) receptor 3B

## Application Notes

Optimal dilution of the Htr3b antibody should be determined by the researcher.

## Immunogen

Amino acids QEDIENDKRSFLND were used as the immunogen for this Htr3b antibody.

## **Storage**

Aliquot and store the Htr3b antibody at -20oC.