


Htr3b Antibody (R35234)

Catalog No.	Formulation	Size
R35234-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	58963
Applications	Western Blot: 0.03-0.1ug/ml ELISA (peptide) LOD: 1:32000
Limitations	This Htr3b antibody is available for research use only.



Description

Additional name(s) for this target protein: 5-hydroxytryptamine (serotonin) receptor 3B

Application Notes

Optimal dilution of the Htr3b antibody should be determined by the researcher.

Immunogen

Amino acids QEDIENDKRSFLND were used as the immunogen for this Htr3b antibody.

Storage Aliquot and store the Htr3b antibody at -20oC.