


# HTR7 Antibody (R35373)

Catalog No.	Formulation	Size
R35373-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3363
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This HTR7 antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: Hydroxytryptamine receptor 7, 5-HT7 receptor

#### **Application Notes**

Optimal dilution of the HTR7 antibody should be determined by the researcher.

1. This HTR7 antibody is specific for isoform a.

# Immunogen

Amino acids QNADYCRKKGHDS were used as the immunogen for this HTR7 antibody.

### **Storage**

Aliquot and store the HTR7 antibody at -20oC.