


HTR7 Antibody (R35373)

Catalog No.	Formulation	Size
R35373-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3363
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This HTR7 antibody is available for research use only.



Western blot testing of KELLY lysate with HTR7 antibody at 1 ug/ml. The expected ~50kDa band and the additional ~37kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Hydroxytryptamine receptor 7, 5-HT7 receptor

Application Notes

Optimal dilution of the HTR7 antibody should be determined by the researcher.

1. This HTR7 antibody is specific for isoform a.

Immunogen

Amino acids QNADYCRKKGHDS were used as the immunogen for this HTR7 antibody.

Storage

Aliquot and store the HTR7 antibody at -20oC.