

Human Immunodeficiency Virus Type-1 p24 Antibody [clone HIV1-24/661] (V2344CF488)

Catalog No.	Formulation	Size
V2344CF488-100T	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 Tests

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	HIV Type 1
Format	CF488 Conjugate
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	HIV1-24/661
Purity	Protein G affinity chromatography
Localization	Membrane
Applications	Immunofluorescence : 1-3ug/ml
Limitations	This Human Immunodeficiency Virus Type-1 p24 antibody is available for research use only.



Description

Human immunodeficiency virus (HIV) is a retrovirus that causes acquired immune deficiency syndrome (AIDS), a condition in humans in which the immune system begins to fail, leading to life-threatening opportunistic infections. HIV mainly infects vital cells in the human immune system such as helper T cells (specifically CD4+ T cells), macrophages and dendritic cells. Two species of HIV infect humans: HIV-1 and HIV-2, with HIV-1 being the more virulent strain. The gag gene of human immunodeficiency virus 1 (HIV-1) encodes a precursor protein known as Pr55Gag. The viral protease PR cleaves this precursor to generate p17, p24, p7, and p6 proteins, which are required for virus particle assembly. HIV-1

Gag p24 is a capsid protein that constitutes the core of AIDS virus HIV-1. p6 and p7 are the components of the nucleocapsid, and p17 provides a protective matrix. HIV-1 Gag p24 is indispensable to the reproduction of AIDS virus and constitutes an essential element for the AIDS virus particle construction. As this protein is detectable from the early stage of AIDS virus infection, its measurement is commonly used as an indicator of HIV-1 infection and viral load.

Application Notes

Optimal dilution of the Human Immunodeficiency Virus Type-1 p24 antibody should be determined by the researcher.

Immunogen

Recombinant HIV-1 Gag p24 protein was used as the immunogen for this Human Immunodeficiency Virus Type-1 p24 antibody.

Storage

Store the Human Immunodeficiency Virus Type-1 p24 antibody at 2-8oC, protected from light.