


IA2 Antibody (R35270)

Catalog No.	Formulation	Size
R35270-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5798
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This IA2 antibody is available for research use only.



Western blot testing of NIH3T3 lysate with IA2 antibody at 0.5ug/ml. Predicted molecular weight ~106kDa.

Description

Additional name(s) for this target protein: Islet cell antigen 512; PTPRN; protein tyrosine phosphatase receptor type N

Application Notes

Optimal dilution of the IA2 antibody should be determined by the researcher.

Immunogen

Amino acids QGEGNIKKNRHPD were used as the immunogen for this IA2 antibody.

Storage

Aliquot and store the IA2 antibody at -20oC.