


ICAM1 Antibody / CD54 (R35091)

Catalog No.	Formulation	Size
R35091-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3383
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ICAM1 antibody is available for research use only.



Western blot testing of HeLa lysate with ICAM1 antibody at 0.1ug/ml. Predicted molecular weight: ~58/75-115kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Intercellular adhesion molecule 1; CD54

Application Notes

Optimal dilution of the ICAM1 antibody should be determined by the researcher.

Immunogen

Amino acids PNVILTKPEVSEGTE were used as the immunogen for this ICAM1 antibody. The amino acid sequence used as immunogen is from the extracellular portion of the protein.

Storage

Aliquot and store the ICAM1 antibody at -20oC.