


ICAM4 Antibody (R35626)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35626-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|--------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 3386 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000 |
| Limitations | This ICAM4 antibody is available for research use only. |



Western blot testing of human umbilical cord lysate with ICAM4 antibody at 1ug/ml.
Predicted molecular weight: ~29kDa.

Description

Additional name(s) for this target protein: Intercellular adhesion molecule 4

Application Notes

Optimal dilution of the ICAM4 antibody should be determined by the researcher.

1. This ICAM4 antibody will detect isoforms 1 and 2.

Immunogen

Amino acids RHGSRVIYSESLER were used as the immunogen for this ICAM4 antibody.

Storage

Aliquot and store the ICAM4 antibody at -20oC.