


IFITM3 Antibody (R36418)

Catalog No.	Formulation	Size
R36418-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10410
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This IFITM3 antibody is available for research use only.



Western blot testing of HeLa cell lysate with IFITM3 antibody at 1ug/ml. Predicted molecular weight: ~15/17kDa (unmodified/lipidated).

Description

Additional name(s) for this target protein: Interferon induced transmembrane protein 3

Application Notes

Optimal dilution of the IFITM3 antibody should be determined by the researcher.

Immunogen

Amino acids HTVQTFFSPVNS were used as the immunogen for this IFITM3 antibody.

Storage

Aliquot and store the IFITM3 antibody at -20oC.