


IFNAR1 Antibody (R35046)

Catalog No.	Formulation	Size
R35046-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3454
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IFNAR1 antibody is available for research use only.



Western blot testing of K562 lysate with IFNAR1 antibody at 0.3ug/ml. Expected molecular weight: 64-135 kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: Interferon alpha receptor 1

Application Notes

Optimal dilution of the IFNAR1 antibody should be determined by the researcher.

Immunogen

Amino acids PKQSGNTPVIQD were used as the immunogen for this IFNAR1 antibody.

Storage

Aliquot and store the IFNAR1 antibody at -20oC.