


IFNAR2 Antibody (RQ5502)


Catalog No.	Formulation	Size
RQ5502	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	P48551
Applications	Western Blot : 0.25-0.5ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml Immunofluorescence : 2-4ug/ml Flow Cytometry : 1-3ug/million cells Direct ELISA : 0.1-0.5ug/ml
Limitations	This IFNAR2 antibody is available for research use only.




Western blot testing of human 1) K562, 2) Caco-2, 3) HepG2 and 4) U-87 MG cell lysate with IFNAR2 antibody. Expected molecular weight: 58-65 kDa.




Immunofluorescent staining of FFPE human A431 with IFNAR2 antibody (green) and DAPI (blue). HIER: boil tissue sections in pH6, 10mM citrate buffer, for 20 min and allow to cool before testing.




IHC staining of FFPE human placenta with IFNAR2 antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 20 min and allow to cool before testing.




IHC staining of FFPE mouse liver with IFNAR2 antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 20 min and allow to cool before testing.



IHC staining of FFPE mouse intestine with IFNAR2 antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 20 min and allow to cool before testing.



IHC staining of FFPE rat intestine with IFNAR2 antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 20 min and allow to cool before testing.



Flow cytometry testing of human A549 cells with IFNAR2 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= IFNAR2 antibody.

Description

Interferon-alpha/beta receptor beta chain is a protein that in humans is encoded by the IFNAR2 gene. It is mapped to 21q22.11. The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinases, which in turn phosphorylate several proteins, including STAT1 and STAT2. Multiple transcript variants encoding at least two different isoforms have been found for this gene.

Application Notes

Optimal dilution of the IFNAR2 antibody should be determined by the researcher.

Immunogen

A human recombinant protein (amino acids D33-Q466) was used as the immunogen for the IFNAR2 antibody.

Storage

After reconstitution, the IFNAR2 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.