


IGFBP7 Antibody (R33812)

Catalog No.	Formulation	Size
R33812-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3490
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This IGFBP7 antibody is available for research use only.



Western blot testing of mouse kidney lysate with IGFBP7 antibody at 1ug/ml. Predicted molecular weight: ~29kDa.

Description

Additional name(s) for this target protein: Insulin-like growth factor binding protein 7

Application Notes

Optimal dilution of the IGFBP7 antibody should be determined by the researcher.

Immunogen

Amino acids EKAITQVSKGTCEQG were used as the immunogen for this IGFBP7 antibody.

Storage

Aliquot and store the IGFBP7 antibody at -20oC.