


IKAROS Antibody (R33536)

Catalog No.	Formulation	Size
R33536-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10320
Localization	Nuclear
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This IKAROS antibody is available for research use only.



Western blot testing of MOLT4 lysate with IKAROS antibody at 0.03ug/ml. Expected molecular weight: ~65/55kDa (isoforms 1/2). (Ref 1)

Description

Additional name(s) for this target protein: IKZF1; IK1

Application Notes

Optimal dilution of the IKAROS antibody should be determined by the researcher.

Immunogen

Amino acids DLLRAASENSQDALR were used as the immunogen for this IKAROS antibody.

Storage

Aliquot and store the IKAROS antibody at -20oC.

References (1)