


IKBIP Antibody (F45794)

Catalog No.	Formulation	Size
F45794-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F45794-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q70UQ0
Applications	Western Blot : 1:1000
Limitations	This IKBIP antibody is available for research use only.



Description

Target of TP53/p53 with pro-apoptotic function.

Application Notes

Titration of the IKBIP antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 319-345 from the human protein was used as the immunogen for this IKBIP antibody.

Storage

Aliquot the IKBIP antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.