


Ikbke Antibody / Ikke / Ikki (F43700)

Catalog No.	Formulation	Size
F43700-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43700-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9R0T8
Applications	Western Blot : 1:1000
Limitations	This Ikbke antibody is available for research use only.



Description

Phosphorylates inhibitors of NF-kappa-B thus leading to the dissociation of the inhibitor/NF-kappa-B complex and ultimately the degradation of the inhibitor. May play a special role in the immune response. Protects cells against DNA damage-induced cell death (By similarity).

Application Notes

Titration of the Ikbke antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 470-498 from the mouse protein was used as the immunogen for this Ikbke antibody. **Storage** Aliquot the Ikbke antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com Copyright $\ensuremath{\texttt{@}}$ NSJ Bioreagents. All rights reserved.