


IL-1 beta Antibody (R35349)

Catalog No.	Formulation	Size
R35349-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3553
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IL-1 beta antibody is available for research use only.



Western blot testing of human peripheral blood lymphocyte lysate with IL-1 beta antibody at 1ug/ml. Predicted molecular weight ~31kDa.

Description

Additional name(s) for this target protein: Interleukin-1 beta; IL1B

Application Notes

Optimal dilution of the IL-1 beta antibody should be determined by the researcher.

Immunogen

Amino acids QLESVDPKNYPKK were used as the immunogen for this IL-1 beta antibody.

Storage

Aliquot and store the IL-1 beta antibody at -20oC.