


IL-17 Antibody (R36188)

Catalog No.	Formulation	Size
R36188-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3605
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This IL-17 antibody is available for research use only.



Western blot testing of human thymus lysate with IL-17 antibody at 1ug/ml. Expected molecular weight: 15-20 kDa (monomer), 30-40 kDa (homodimer), depending on glycosylation level.

Description

Additional name(s) for this target protein: Interleukin 17A antibody, IL17A, CTLA8

Application Notes

Optimal dilution of the IL-17 antibody should be determined by the researcher.

Immunogen

Amino acids TIPRNPGCPNSEDKN from human IL-17A were used as the immunogen for this IL-17 antibody.

Storage

Aliquot and store the IL-17 antibody at -20oC.