

IL-18 Antibody (R34876)

Catalog No.	Formulation	Size
R34876-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3606
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This IL-18 antibody is available for research use only.



Western blot testing of HeLa lysate with IL-18 antibody at 0.1ug/ml. Predicted molecular weight ~23kDa.

Description

Additional name(s) for this target protein: Interleukin 18, interleukin-1 gamma; IL18

Application Notes

Optimal dilution of the IL-18 antibody should be determined by the researcher.

Immunogen

Amino acids NPPDNIKDTKSDI were used as the immunogen for this IL-18 antibody.

Storage

Aliquot and store the IL-18 antibody at -20oC.