


IL-20 Antibody / IL-17B (R35879)

Catalog No.	Formulation	Size
R35879-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27190
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IL-20 antibody is available for research use only.



Western blot testing of human Ileum lysate with IL-20 antibody at 2ug/ml. The ~20 kDa (glycosylated) band and the ~17 kDa (unmodified) band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Interleukin 20, interleukin 17B, IL17B; IL17B

Application Notes

Optimal dilution of the IL-20 antibody should be determined by the researcher.

Immunogen

Amino acids SELAQRKCEVNLQ were used as the immunogen for this IL-20 antibody.

Storage

Aliquot and store the IL-20 antibody at -20oC.