


IL-6 Antibody (R31438)

Catalog No.	Formulation	Size
R31438	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide
UniProt	P20607
Applications	Western Blot : 0.5-1ug/ml
Limitations	This IL-6 antibody is available for research use only.



Western blot testing of IL-6 antibody (0.5ug/ml) and recombinant rat protein (0.5ng)

Description

Interleukin-6 is a protein that in humans is encoded by the IL6 gene. It is an interleukin that acts as both a pro-inflammatory and anti-inflammatory cytokine. It is secreted by T cells and macrophages to stimulate immune response to trauma, especially burns or other tissue damage leading to inflammation. IL-6 is one of the most important mediators of fever and of the acute phase response. It is also essential for hybridoma growth and is found in many supplemental cloning media such as briclone. Bowcock et al. (1988) assigned the IL6 gene to chromosome 7p21. By in situ hybridization and Southern blot analysis of mouse-human hybrid cell lines, Sutherland et al. (1988) mapped the gene to chromosome 7p15.

Application Notes

The stated application concentrations are suggested starting points. Titration of the IL-6 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Rat partial recombinant protein (AA 25-211) was used as the immunogen for this IL-6 antibody.

Storage

After reconstitution, the IL-6 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.