


## Il12b Antibody (R33713)

Catalog No.	Formulation	Size
R33713-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	16160
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Il12b antibody is available for research use only.



Western blot testing of 1) mouse skin and 2) rat skin lysate with Il12b antibody at 0.5ug/ml. Predicted molecular weight: 37~40 kDa.

## Description

Additional name(s) for this target protein: Interleukin 12B

## Application Notes

Optimal dilution of the Il12b antibody should be determined by the researcher.

## Immunogen

Amino acids EARQQNKYENYS were used as the immunogen for this IL12b antibody.

## **Storage**

Aliquot and store the IL12b antibody at -20°C.