


## IL12B Antibody (R35275)

Catalog No.	Formulation	Size
R35275-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3593
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IL12B antibody is available for research use only.



Western blot testing of human 1) liver, 2) skin and 3) tonsil lysate with IL12B antibody at 0.3ug/ml. Predicted molecular weight: 37~40 kDa.

## Description

Additional name(s) for this target protein: Interleukin 12B (natural killer cell stimulatory factor 2, cytotoxic lymphocyte maturation factor 2, p40); IL12/p40

## Application Notes

Optimal dilution of the IL12B antibody should be determined by the researcher.

## **Immunogen**

Amino acids QGKSKREKKDRVFTD were used as the immunogen for this IL12B antibody.

## **Storage**

Aliquot and store the IL12B antibody at -20oC.