

## IL12RB1 Antibody / CD212 (R34322)

Catalog No.	Formulation	Size
R34322-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3594
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This IL12RB1 antibody is available for research use only.



Western blot testing of peripheral blood lymphocyte lysate with IL12RB1 antibody at 0.3ug/ml. Predicted molecular weight: ~73/100kDa (unmodified/glycosylated).

## Description

Additional name(s) for this target protein: Interleukin 12 receptor, beta 1; CD212

## Application Notes

Optimal dilution of the IL12RB1 antibody should be determined by the researcher.

1. This IL12RB1 antibody is specific for isoform 1.

## Immunogen

Amino acids HVSVKNHSLDSVS were used as the immunogen for this IL12RB1 antibody.

## Storage

Aliquot and store the IL12RB1 antibody at -20oC.