


IL12RB1 Antibody / CD212 (R34322)

Catalog No.	Formulation	Size
R34322-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3594
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This IL12RB1 antibody is available for research use only.



Western blot testing of peripheral blood lymphocyte lysate with IL12RB1 antibody at 0.3ug/ml. Predicted molecular weight: ~73/100kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Interleukin 12 receptor, beta 1; CD212

Application Notes

Optimal dilution of the IL12RB1 antibody should be determined by the researcher.

1. This IL12RB1 antibody is specific for isoform 1.

Immunogen

Amino acids HVSVKNHSLDSVS were used as the immunogen for this IL12RB1 antibody.

Storage

Aliquot and store the IL12RB1 antibody at -20oC.