


## IL13RA2 Antibody (R33538)

Catalog No.	Formulation	Size
R33538-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3598
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IL13RA2 antibody is available for research use only.



Western blot testing of HeLa lysate with IL13RA2 antibody at 1ug/ml. Predicted molecular weight: ~44/60kDa (unmodified/glycosylated).

## Description

Additional name(s) for this target protein: Interleukin 13 receptor alpha 2

## Application Notes

Optimal dilution of the IL13RA2 antibody should be determined by the researcher.

## Immunogen

Amino acids DLNKGIEAKIHT were used as the immunogen for this IL13RA2 antibody.

## Storage

Aliquot and store the IL13RA2 antibody at -20oC.

## References (1)