


IL15 Antibody (R33386)

Catalog No.	Formulation	Size
R33386-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3600
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IL15 antibody is available for research use only.



Western blot testing of human bone marrow lysate with IL15 antibody at 1ug/ml.
Predicted molecular weight: ~18/23kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Interleukin 15

Application Notes

Optimal dilution of the IL15 antibody should be determined by the researcher.

Immunogen

Amino acids DATLYTESDVHPS were used as the immunogen for this IL15 antibody.

Storage

Aliquot and store the IL15 antibody at -20oC.