


IL17C Antibody (R35752)

Catalog No.	Formulation	Size
R35752-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27189
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IL17C antibody is available for research use only.



Western blot testing of human lymph node lysate with IL17C antibody at 0.3ug/ml.
Predicted molecular weight: ~22 kDa.

Description

Additional name(s) for this target protein: Interleukin 17C

Application Notes

Optimal dilution of the IL17C antibody should be determined by the researcher.

Immunogen

Amino acids RRPCSRDGSGGLPT were used as the immunogen for this IL17C antibody.

Storage

Aliquot and store the IL17C antibody at -20oC.