


IL17D Antibody (F44066)

Catalog No.	Formulation	Size
F44066-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F44066-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8TAD2
Applications	Western Blot : 1:1000
Limitations	This IL17D antibody is available for research use only.



Description

The protein encoded by this gene is a cytokine that shares the sequence similarity with IL17. The treatment of endothelial cells with this cytokine has been shown to stimulate the production of other cytokines including IL6, IL8 and CSF2/ GM-CSF. The increased expression of IL8 induced by this cytokine was found to be NF-kappa B-dependent.

Application Notes

Titration of the IL17D antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 155-183 from the human protein was used as the immunogen for this IL17D antibody. **Storage** Aliquot the IL17D antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com Copyright $\ensuremath{\texttt{@}}$ NSJ Bioreagents. All rights reserved.