


IL17ra Antibody (R33591)

Catalog No.	Formulation	Size
R33591-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	16172
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This IL17ra antibody is available for research use only.



Western blot testing of rat kidney lysate with IL17ra antibody at 1ug/ml. Predicted molecular weight ~97 kDa but can be observed at a larger size due to glycosylation.

Description

Additional name(s) for this target protein: Interleukin 17 receptor A

Application Notes

Optimal dilution of the IL17ra antibody should be determined by the researcher.

Immunogen

Amino acids DQEKHGDDSKIN were used as the immunogen for this IL17ra antibody.

Storage

Aliquot and store the IL17ra antibody at -20oC.