


IL17RD Antibody (F51826)

| Catalog No. | Formulation | Size |
|---------------|--|---------|
| F51826-0.4ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.4 ml |
| F51826-0.08ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.08 ml |

[Bulk quote request](#)

| | |
|---------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (rabbit origin) |
| Isotype | Rabbit Ig |
| Purity | Antigen affinity |
| UniProt | Q8NFM7 |
| Applications | Western Blot : 1:1000 |
| Limitations | This IL17RD antibody is available for research use only. |



Description

Fibroblast growth factors are secreted proteins involved in cellular proliferation, migration, differentiation, and survival. FGF activity is negatively regulated by members of the 'sprouty' family. SEF/IL17RD is a modulator of FGF signaling.

Application Notes

Titration of the IL17RD antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 217-245 from the human protein was used as the immunogen for this IL17RD antibody.

Storage

Aliquot the IL17RD antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.