

IL33 Antibody (R34000)

Catalog No.	Formulation	Size
R34000-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	90865
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This IL33 antibody is available for research use only.



Western blot testing of Jurkat lysate with IL33 antibody at 1ug/ml. Predicted molecular weight: ~31/26/16kDa (isoforms a/b/c).

Description

Additional name(s) for this target protein: Interleukin 33; IL-33

Application Notes

Optimal dilution of the IL33 antibody should be determined by the researcher.

Immunogen

Amino acids ESQHPSNESGDGVD were used as the immunogen for this IL33 antibody. This sequence is common to isoforms a/b/c.

Storage

Aliquot and store the IL33 antibody at -20oC.