


## Importin 7 Antibody / IPO7 (R36066)

Catalog No.	Formulation	Size
R36066-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Mouse, Rat, Dog, Cow
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	10527
<b>Applications</b>	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:32000
<b>Limitations</b>	This Importin 7 antibody is available for research use only.



Western blot testing of Daudi lysate with Importin 7 antibody at 1 ug/ml. The expected ~120 kDa band and the additional ~26kDa band are both blocked by the immunizing peptide.

### Description

Additional name(s) for this target protein: IPO7, RAN binding protein 7, RANBP7

### Application Notes

Optimal dilution of the Importin 7 antibody should be determined by the researcher.

## **Immunogen**

Amino acids SSFNFGGPAPGMN were used as the immunogen for this Importin 7 antibody.

## **Storage**

Aliquot and store the Importin 7 antibody at -20oC.