


## Importin 7 Antibody / IPO7 (R36066)

Catalog No.	Formulation	Size
R36066-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10527
Applications	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Importin 7 antibody is available for research use only.



Western blot testing of Daudi lysate with Importin 7 antibody at 1 ug/ml. The expected ~120 kDa band and the additional ~26kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: IPO7, RAN binding protein 7, RANBP7

## Application Notes

Optimal dilution of the Importin 7 antibody should be determined by the researcher.

## Immunogen

Amino acids SSFNFGGPAPGMN were used as the immunogen for this Importin 7 antibody.

## **Storage**

Aliquot and store the Importin 7 antibody at -20oC.