


INCA Antibody (R34818)

Catalog No.	Formulation	Size
R34818-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	440068
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This INCA antibody is available for research use only.



Western blot testing of human kidney lysate with INCA antibody at 0.5ug/ml. Predicted molecular weight: ~12kDa.

Description

Additional name(s) for this target protein: Caspase recruitment domain family, member 17, CARD17, Caspase-1 inhibitor INCA

Application Notes

Optimal dilution of the INCA antibody should be determined by the researcher.

Immunogen

Amino acids NHLTTQDSQIVLPS were used as the immunogen for this INCA antibody.

Storage

Aliquot and store the INCA antibody at -20oC.