


ING2 Antibody / Inhibitor of growth protein 2 (R33513)

Catalog No.	Formulation	Size
R33513-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3622
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This ING2 antibody is available for research use only.



Western blot testing of HepG2 lysate with ING2 antibody at 2ug/ml. Predicted molecular weight: ~33kDa.

Description

Additional name(s) for this target protein: Inhibitor of growth family, member 2, p33ING2

Application Notes

Optimal dilution of the ING2 antibody should be determined by the researcher.

Immunogen

Amino acids DKSTKTKKDRRSR were used as the immunogen for this ING2 antibody.

Storage

Aliquot and store the ING2 antibody at -20oC.