


Integrin alpha 1 Antibody (R33502)

Catalog No.	Formulation	Size
R33502-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3672
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Integrin alpha 1 antibody is available for research use only.



Western blot testing of HeLa [A] and HepG2 [B] lysates with Integrin alpha 1 antibody at 0.5ug/ml. Predicted molecular weight: ~131/150+ kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: ITGA1; CD49; CD49a, VLA1

Application Notes

Optimal dilution of the Integrin alpha 1 antibody should be determined by the researcher.

Immunogen

Amino acids DKHDFQDSVRIT were used as the immunogen for this Integrin alpha 1 antibody.

Storage

Aliquot and store the Integrin alpha 1 antibody at -20oC.