


IP6K3 Antibody / IHPK3 (R34419)

Catalog No.	Formulation	Size
R34419-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	117283
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IP6K3 antibody is available for research use only.



Western blot testing of human kidney lysate with IP6K3 antibody at 1ug/ml.
Predicted/observed molecular weight: ~46kDa.

Description

Additional name(s) for this target protein: Inositol hexakisphosphate kinase 3

Application Notes

Optimal dilution of the IP6K3 antibody should be determined by the researcher.

Immunogen

Amino acids KESQEPFKVSTES were used as the immunogen for this IP6K3 antibody.

Storage

Aliquot and store the IP6K3 antibody at -20oC.