


IQGAP1 Antibody (R34633)

Catalog No.	Formulation	Size
R34633-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8826
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This IQGAP1 antibody is available for research use only.



Western blot testing of human kidney lysate with IQGAP1 antibody at 1ug/ml. Predicted molecular weight ~189 kDa.

Description

Additional name(s) for this target protein: IQ motif containing GTPase activating protein 1

Application Notes

Optimal dilution of the IQGAP1 antibody should be determined by the researcher.

Immunogen

Amino acids KRQSGQTDPLQKEE were used as the immunogen for this IQGAP1 antibody.

Storage

Aliquot and store the IQGAP1 antibody at -20oC.