


Irak2 Antibody (F44544)

Catalog No.	Formulation	Size
F44544-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F44544-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8CFA1
Applications	Western Blot : 1:1000
Limitations	This Irak2 antibody is available for research use only.



Irak2 antibody western blot analysis in A375 lysate.

Description

Irak2 binds to the IL-1 type I receptor following IL-1 engagement, triggering intracellular signaling cascades leading to transcriptional up-regulation and mRNA stabilization (By similarity).

Application Notes

Titration of the Irak2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 594-622 from the mouse protein was used as the immunogen for this Irak2 antibody.

Storage

Aliquot the Irak2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.