


Irak3 Antibody (F43676)

Catalog No.	Formulation	Size
F43676-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43676-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8K4B2
Applications	Western Blot : 1:1000
Limitations	This Irak3 antibody is available for research use only.



Description

Irak3 inhibits dissociation of IRAK1 and IRAK4 from the Toll-like receptor signaling complex by either inhibiting the phosphorylation of IRAK1 and IRAK4 or stabilizing the receptor complex.

Application Notes

Titration of the Irak3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 517-544 from the mouse protein was used as the immunogen for this Irak3 antibody.

Storage

Aliquot the Irak3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.