

## IRAK3 Antibody (R34098)


Catalog No.	Formulation	Size
R34098-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	73914
Applications	Western Blot : 0.05-0.2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This IRAK3 antibody is available for research use only.



Western blot testing of rat liver lysate with IRAK3 antibody at 0.05ug/ml. Predicted molecular weight: ~68kDa.



Western blot of mouse bone marrow (wildtype first lane, KO second lane) lysates with IRAK3 antibody at 0.5ug/ml. Predicted molecular weight: ~68kDa.

## Description

Additional name(s) for this target protein: Interleukin-1 receptor-associated kinase 3, Interleukin-1 receptor-associated kinase M

## Application Notes

Optimal dilution of the IRAK3 antibody should be determined by the researcher.

## Immunogen

Amino acids GHSYGSKPMEKR were used as the immunogen for this IRAK3 antibody.

## Storage

Aliquot and store the IRAK3 antibody at -20oC.