


IRF4 Antibody (R34209)


Catalog No.	Formulation	Size
R34209-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3662
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This IRF4 antibody is available for research use only.



Western blot testing of Jurkat lysate with IRF4 antibody at 0.5ug/ml. Predicted molecular weight: 51 kDa.



Western blot testing of human 1) Daudi and 2) Jurkat lysate with IRF4 antibody at 0.5ug/ml. Predicted molecular weight: 51 kDa.

Description

Additional name(s) for this target protein: Interferon regulatory factor 4

Application Notes

Optimal dilution of the IRF4 antibody should be determined by the researcher.

Immunogen

Amino acids HISNPEDYHRSIR were used as the immunogen for this IRF4 antibody.

Storage

Aliquot and store the IRF4 antibody at -20oC.